B. Braun Medical AG

Meliseptol Foam pure / fresh

Revision date: 01.06.2015 Revision No: 1,1

Product code: 00056-0328

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1. Product identifier

Meliseptol Foam pure / fresh

#### 1.2. Relevant identified uses of the substance or mixture and uses advised against

#### Use of the substance/mixture

Rapid disinfectant for small surfaces

### 1.3. Details of the supplier of the safety data sheet

Manufacturer

Company name: B. Braun Medical AG

Seesatz 17 Street:

Place: CH-6204 Sempach

Responsible Department: Telephone: +41 (0) 58 / 258 50 00

E-Mail: info.bbmch@bbraun.com Responsible for the safety data sheet: sds@gbk-ingelheim.de

**Supplier** 

Company name: B. Braun Melsungen AG Carl-Braun-Straße 1 Street: Place: D-34212 Melsungen

Responsible Department: Zentrale Service-Bereiche / Logistik und Supply Chain

> Telephone: +49 (0) 5661 / 71-4422 E-Mail: logistics.service@bbraun.com

1.4. Emergency telephone number: INTERNATIONAL: +49 - (0) 6132 - 84463, GBK GmbH (24h - 7d/w - 365d/a) England and Wales: NHS Direct - 0845 4647; Scotland: NHS 24 - 08454 24 24 24

#### SECTION 2: Hazards identification

#### 2.1. Classification of the substance or mixture according to 1272/2008/EC

Hazard categories:

Flammable liquid: Flam. Liq. 3

Serious eye damage/eye irritation: Eye Dam. 1

Hazard Statements:

Flammable liquid and vapour. Causes serious eye damage.

#### 2.2. Label elements

## Hazardous components which must be listed on the label

Propan-1-ol

Signal word: Danger

Pictograms:



#### **Hazard statements**

H226 Flammable liquid and vapour. H318 Causes serious eye damage.

## **Precautionary statements**

Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No P210

smoking.

P260 Do not breathe vapour.

Wear protective gloves/protective clothing/eye protection/face protection. P280

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if

present and easy to do. Continue rinsing.

IF INHALED: Remove person to fresh air and keep comfortable for breathing. P304+P340

Call a POISON CENTER/doctor if you feel unwell. P312

> GB - FN Page 1 of 8

B. Braun Medical AG

## Meliseptol Foam pure / fresh

Revision date: 01.06.2015 Revision No: 1,1

Product code: 00056-0328

#### 2.3. Other hazards

In use, may form flammable/explosive vapour-air mixture.

## **SECTION 3: Composition/information on ingredients**

#### 3.2. Mixtures

#### **Chemical characterization**

Alcoholic solution

#### **Hazardous components**

CAS No	Chemical name				
	EC No	Index No	REACH No		
	Classification according to Regulation (EC) No. 1272/2008 [CLP]				
71-23-8	Propan-1-ol			< 20 %	
	200-746-9	603-003-00-0	01-2119486761-29		
	Flam. Liq. 2, Eye Dam. 1,				
7173-51-5	Didecyldimethylammonium chloride			< 0,25 %	
	230-525-2	612-131-00-6	01-2119945987-15		
	Acute Tox. 4, Skin Corr. 1B, Aquatic Acute 1 (M-Factor = 10); H302 H314 H400				

Full text of H and EUH statements: see section 16.

## **SECTION 4: First aid measures**

## 4.1. Description of first aid measures

#### **General information**

Remove contaminated soaked clothing immediately.

If you feel unwell, seek medical advice.

#### After inhalation

Move to fresh air in case of accidental inhalation of vapours.

In the event of symptoms refer for medical treatment.

### After contact with skin

In case of contact with skin wash off immediately with plenty of water.

Consult a doctor if skin irritation persists.

#### After contact with eyes

Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes.

Seek medical treatment by eye specialist.

## After ingestion

Drink plenty of water.

Summon a doctor immediately.

Induce vomiting only upon the advice of a physician.

## 4.2. Most important symptoms and effects, both acute and delayed

Causes serious eye damage.

## 4.3. Indication of any immediate medical attention and special treatment needed

Treat symptoms.

## **SECTION 5: Firefighting measures**

## 5.1. Extinguishing media

## Suitable extinguishing media

Alcohol-resistant foam, dry chemical, carbon dioxide (CO2), water-spray.

## Unsuitable extinguishing media

Full water jet.

### 5.2. Special hazards arising from the substance or mixture

Fire may produce:

Carbon monoxide (CO), carbon dioxide (CO2) and nitrogen oxides (NOx).

Chlorine compounds.

GB - EN Page 2 of 8

B. Braun Medical AG

#### Meliseptol Foam pure / fresh

Revision date: 01.06.2015 Revision No: 1,1

Product code: 00056-0328

## 5.3. Advice for firefighters

Use breathing apparatus with independent air supply.

Protective suit.

#### Additional information

Cool containers at risk with water spray jet.

Keep away from sources of ignition - No smoking.

Vapours are heavier than air and spread along ground.

The vapour/air mixture is explosive, even in empty, uncleaned receptacles.

## **SECTION 6: Accidental release measures**

## 6.1. Personal precautions, protective equipment and emergency procedures

In case of vapour formation use respirator.

Avoid contact with skin, eyes and clothing.

Ensure adequate ventilation.

Keep away sources of ignition.

#### 6.2. Environmental precautions

Do not discharge into the drains/surface waters/ground water.

## 6.3. Methods and material for containment and cleaning up

Soak up with inert absorbent material (e.g. sand, silica gel, acid binder, universal binder).

Shovel into suitable container for disposal.

#### 6.4. Reference to other sections

Observe protective instructions (see Sections 7 and 8).

Information for disposal see section 13.

## **SECTION 7: Handling and storage**

## 7.1. Precautions for safe handling

### Advice on safe handling

Ensure adequate ventilation.

When using do not eat, drink or smoke.

Avoid contact with eyes.

### Advice on protection against fire and explosion

Keep product and empty container away from heat and sources of ignition.

Do not smoke - volatile.

Take precautionary measures against static discharges.

## 7.2. Conditions for safe storage, including any incompatibilities

### Requirements for storage rooms and vessels

Keep container tightly closed in a dry, cool and well-ventilated place.

Pay attention to anti-explosion rules.

#### Advice on storage compatibility

Incompatible with:

Oxidizing agents

Alkaline metals and earth alkaline metals.

#### Further information on storage conditions

Keep away from food, drink and animal feeding stuffs.

#### 7.3. Specific end use(s)

Rapid disinfectant for small surfaces

## **SECTION 8: Exposure controls/personal protection**

## 8.1. Control parameters

## **Exposure limits (EH40)**

CAS No	Substance	ppm	mg/m³	fibres/ml	Category	Origin
71-23-8	Propan-1-ol	200	500		TWA (8 h)	WEL
		250	625		STEL (15 min)	WEL

B. Braun Medical AG

Meliseptol Foam pure / fresh

Revision date: 01.06.2015 Revision No: 1,1

Product code: 00056-0328

#### 8.2. Exposure controls

### Appropriate engineering controls

Ensure adequate ventilation, especially in confined areas.

## Protective and hygiene measures

Do not inhale vapours.

When using do not eat, drink or smoke.

Avoid contact with the eyes.

Remove and wash contaminated clothing before re-use.

#### Eye/face protection

Safety goggles with side protection (EN 166).

Eye wash bottle with pure water (EN 15154).

#### Hand protection

Protective gloves resistant to chemicals made of nitrile, minimum coat thickness 0.4 mm, permeation resistance (wear duration) approx. 480 minutes, i.e. protective glove <Camatril Velours 730> made by www.kcl.de.

This recommendation refers exclusively to the chemical compatibility and the lab test conforming to EN 374 carried out under lab conditions.

Requirements can vary as a function of the use. Therefore it is necessary to adhere additionally to the recommendations given by the manufacturer of protective gloves.

#### Respiratory protection

No personal respiratory protective equipment normally required.

In case of insufficient ventilation wear suitable respiratory equipment (gas filter type A) (EN 14387).

## **SECTION 9: Physical and chemical properties**

#### 9.1. Information on basic physical and chemical properties

Physical state: Liquid
Colour: Colourless
Odour: characteristic

pH-Value (at 20 °C): approx. 7

Changes in the physical state

Initial boiling point and boiling range: approx. 89 °C

Flash point: 31,5 °C DIN 51755 \*)

Lower explosion limits: n.d. Ignition temperature: n.d.

Density (at 20 °C): approx. 0,975 g/cm<sup>3</sup>

Water solubility: Miscible

(at 20 °C)

Solvent content: < 20 %

#### 9.2. Other information

\*) Product is ignitable, but does not keep burning.

## **SECTION 10: Stability and reactivity**

## 10.1. Reactivity

No decomposition if stored and applied as directed.

## 10.2. Chemical stability

Stable under normal conditions.

## 10.3. Possibility of hazardous reactions

Reactions with oxidizing agents.

Reactions with alkali metals.

Reactions with earth alkali metals.

#### 10.4. Conditions to avoid

Vapour/air mixtures are explosive at intensive warming.

GB - EN Page 4 of 8

B. Braun Medical AG

Meliseptol Foam pure / fresh

Revision date: 01.06.2015 Revision No: 1,1

Product code: 00056-0328

Heating can release vapours which can be ignited.

#### 10.5. Incompatible materials

oxidizing agents

Alkaline metals and alkaline earth metals.

#### 10.6. Hazardous decomposition products

Carbon monoxide (CO), carbon dioxide (CO2) and nitrogen oxides (NOx).

Chlorine compounds.

## **SECTION 11: Toxicological information**

#### 11.1. Information on toxicological effects

#### **Acute toxicity**

Based on available data, the classification criteria are not met.

No toxicological data available.

## Irritation and corrosivity

Causes serious eye damage. Skin irritation: Not classified.

#### Sensitising effects

Based on available data, the classification criteria are not met.

#### STOT-single exposure

Based on available data, the classification criteria are not met.

#### Severe effects after repeated or prolonged exposure

Based on available data, the classification criteria are not met.

### Carcinogenic/mutagenic/toxic effects for reproduction

Based on available data, the classification criteria are not met.

#### **Aspiration hazard**

Based on available data, the classification criteria are not met.

#### Additional information on tests

Classification in compliance with the assessment procedure specified in the Regulation (EC) no 1272/2008.

## **Practical experience**

#### Other observations

Ingestion causes irritation of upper respiratory system and gastrointestinal disturbance. Swallowing renders reabsorption possible.

Repeated or prolonged exposure may cause skin irritation and dermatitis, due to degreasing properties of the product.

Inhalation of vapours in high concentration can cause narcotic effects.

## **SECTION 12: Ecological information**

## 12.1. Toxicity

Ecological data are not available.

Didecyldimethylammonium chloride [M = 10]

LC50/Pimephales promelas/96 h = 0,19 mg/l [US-EPA]

ErC50/Pseudokirchneriela subcapitata/96 h = 0,026 mg/l [OECD TG 201]

EC50/Daphnia magna/48 h = 0,062 mg/l [EPA-FIFRA]

#### 12.2. Persistence and degradability

The surfactant(s) contained in this preparation complies(comply) with the biodegradability criteria as laid down in Regulation (EC) No.648/2004 on detergents. Data to support this assertion are held at the disposal of the competent authorities of the Member States and will be made available to them, at their direct request or at the request of a detergent manufacturer.

CAS No	Chemical name				
	Method	Value	d	Source	
	Evaluation				
7173-51-5	Didecyldimethylammonium chloride				
	modif. Sturm-test	72 %	28		
	Readily biodegradable.				

B. Braun Medical AG

Meliseptol Foam pure / fresh

Revision date: 01.06.2015 Revision No: 1,1

Product code: 00056-0328

## 12.3. Bioaccumulative potential

Propan-1-ol

Product has a low bioaccumulating potential.

### 12.4. Mobility in soil

No data available.

#### 12.5. Results of PBT and vPvB assessment

According to Regulation (EC) No 1907/2006 (REACH) none of the substances, contained in this product are a PBT substance.

#### 12.6. Other adverse effects

Low hazard to waters.

#### **Further information**

Do not flush into surface water or sanitary sewer system.

When low concentrations are discharged correctly into adapted biological sewage treatment plants. interference with the degradation activity of activated sludge is not likely.

## **SECTION 13: Disposal considerations**

### 13.1. Waste treatment methods

#### Advice on disposal

Can be disposed of as a solid waste or burned in a suitable installation subject to local regulations. Where possible recycling is preferred to disposal.

## Waste disposal number of waste from residues/unused products

070604 WASTES FROM ORGANIC CHEMICAL PROCESSES; wastes from the MFSU of fats, grease, soaps,

detergents, disinfectants and cosmetics; other organic solvents, washing liquids and mother liquors

Classified as hazardous waste.

#### Contaminated packaging

Empty containers should be taken for local recycling, recovery or waste disposal.

Contaminated packaging should be emptied as far as possible and after appropriate cleansing may be taken for reuse.

Packaging that cannot be cleaned should be disposed of like the product.

## **SECTION 14: Transport information**

## Land transport (ADR/RID)

14.1. UN number: UN 1274

14.2. UN proper shipping name: n-PROPANOL, SOLUTION

3 14.3. Transport hazard class(es): 14.4. Packing group: Ш 3

Hazard label:



Classification code:

5 L / 30 kg Limited quantity:

E1 Excepted quantity: Transport category: 3 Hazard No: 30 Tunnel restriction code: D/E

Inland waterways transport (ADN)

14.1. UN number: UN 1274

14.2. UN proper shipping name: n-PROPANOL, SOLUTION

14.3. Transport hazard class(es): 3 Ш 14.4. Packing group: Hazard label: 3

> GB - FN Page 6 of 8

B. Braun Medical AG

Meliseptol Foam pure / fresh

Revision date: 01.06.2015 Revision No: 1,1

Product code: 00056-0328



Classification code: F1

Limited quantity: 5 L / 30 kg

Excepted quantity: E1

Marine transport (IMDG)

**14.1. UN number:** UN 1274

**14.2. UN proper shipping name:** n-PROPANOL (PROPYL ALCOHOL, NORMAL), SOLUTION

14.3. Transport hazard class(es):314.4. Packing group:IIIHazard label:3



Marine pollutant: No

Limited quantity: 5 L / 30 kg

Excepted quantity: E1

EmS: F-E, S-D

Air transport (ICAO)

**14.1. UN number:** UN 1274

**14.2. UN proper shipping name:** n-PROPANOL (PROPYL ALCOHOL, NORMAL), SOLUTION

14.3. Transport hazard class(es):314.4. Packing group:IIIHazard label:3



Limited quantity Passenger: 10 L
Passenger LQ: Y344
Excepted quantity: E1

IATA-packing instructions - Passenger:355IATA-max. quantity - Passenger:60 LIATA-packing instructions - Cargo:366IATA-max. quantity - Cargo:220 L

14.5. Environmental hazards

ENVIRONMENTALLY HAZARDOUS: no

14.6. Special precautions for user

Handle in accordance with good industrial hygiene and safety practice.

14.7. Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code

The transport takes place only in approved and appropriate packaging.

# **SECTION 15: Regulatory information**

## 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

EU regulatory information

2004/42/EC (VOC): < 20 %

Additional information

Regulation (EC) No 648/2004 (Regulation on detergents):

Non-ionic surfactants < 5 %, Perfume < 5% (Meliseptol Foam fresh)

Ingredients subject to the labelling obligation according to SCCP: -

## National regulatory information

GB - EN Page 7 of 8

B. Braun Medical AG

Meliseptol Foam pure / fresh

Revision date: 01.06.2015 Revision No: 1,1

Product code: 00056-0328

Employment restrictions: Observe employment restrictions for young people. Observe employment

restrictions for child bearing mothers and nursing.

Water contaminating class (D): 1 - slightly water contaminating

#### 15.2. Chemical safety assessment

For this substance a chemical safety assessment has not been carried out.

#### **SECTION 16: Other information**

#### Abbreviations and acronyms

ADR = Accord européen relatif au transport international des marchandises Dangereuses par Route

RID = Règlement concernant le transport international ferroviaire de marchandises dangereuses

ADN = Accord européen relatif au transport international des marchandises dangereuses par voie de navigation intérieure IMDG = International Maritime Code for Dangerous Goods

IATA/ICAO = International Air Transport Association / International Civil Aviation Organization

MARPOL = International Convention for the Prevention of Pollution from Ships

IBC-Code = International Code for the Construction and Equipment of Ships Carrying Dangerous Chemicals in Bulk

GHS = Globally Harmonized System of Classification and Labelling of Chemicals

REACH = Registration, Evaluation, Authorization and Restriction of Chemicals

CAS = Chemical Abstract Service

EN = European norm

ISO = International Organization for Standardization

DIN = Deutsche Industrie Norm

PBT = Persistent Bioaccumulative and Toxic

vPvB = Very Persistent and very Bio-accumulative

LD = Lethal dose

LC = Lethal concentration

EC = Effect concentration

IC = Median immobilisation concentration or median inhibitory concentration

#### Relevant H and EUH statements (number and full text)

H225	Highly flammable liquid and vapour.
H226	Flammable liquid and vapour.
H302	Harmful if swallowed.
H314	Causes severe skin burns and eye damage.
H318	Causes serious eve damage

H318 Causes serious eye damage.
H336 May cause drowsiness or dizziness.

H400 Very toxic to aquatic life.

#### **Further Information**

Data of items 4 to 8, as well as 10 to 12, do partly not refer to the use and the regular employing of the product (in this sense consult information on use and on product), but to liberation of major amounts in case of accidents and irregularities.

The information describes exclusively the safety requirements for the product(s) and is based on the present level of our knowledge.

The delivery specifications are contained in the corresponding product sheet.

This data does not constitute a guarantee for the characteristics of the product(s) as defined by the legal warranty regulations.

(n.a. = not applicable; n.d. = not determined)

(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)

GB - EN Page 8 of 8